

α 2-Antiplasmin Activity Assay

Antiplasmin activity is measured by level of inhibition of plasmin

Materials

1. Versa Max microtiter plate reader (Molecular Devices, CA)
2. 96-well microtiter plate
3. Citrated plasma samples
4. Normal mouse plasma (NMP)
5. TBS with 120 mM methylamine chloride
6. 10 μ g/mL plasmin
7. 1.5 mM chromogenic substrate specific for plasmin
8. 20% acetic acid

Procedure

1. Plasma samples may be diluted to 1:50 in TBS/120 mM methylamine chloride
2. Standard curve is prepared with each plate by diluting NMP 1:10 to 1:1280 in TBS/120 mM methylamine chloride, analyzed simultaneously on the sample plate
3. Plate wells are loaded as follows:

Sample dilutions, in duplicate	50uL	37°C	15min	
Plasmin	50uL		90sec	
Chromogenic substrate	50uL		10min	Color development
Acetic acid	50uL			Stop reaction

4. Plate is read at 405 nm at desired time points as color develops and after the reaction is stopped
5. Absorbance are converted to %NMP α 2-antiplasmin using the standard curve