

## Plasminogen Activity Assay

Activity is measured by the level of chromogenic substrate cleaved by plasmin or urokinase-activated plasminogen

### Materials

1. Versa Max microtiter plate reader (Molecular Devices, CA)
2. 96-well microtiter plate
3. Citrated plasma samples
4. Normal mouse plasma (NMP)
5. 100 mM Tris, pH 8.5 with 8.3 mM EACA (Tris/EACA)
6. HN/BSA
7. 2500 Ploug U/mL urokinase
8. 1.2 mM chromogenic substrate specific for plasmin and urokinase-activated plasminogen in HN buffer
9. 20% acetic acid

### Procedure

1. Plasma samples may be diluted to 1:50 in Tris/EACA
2. Log-log standard curve is prepared with each plate by diluting NMP 1:10 to 1:320 in Tris/EACA, analyzed simultaneously on the sample plate
3. Plate wells are loaded as follows:

|                                       |              |             |              |                          |
|---------------------------------------|--------------|-------------|--------------|--------------------------|
| <b>Sample dilutions, in duplicate</b> | <b>60uL</b>  | <b>37°C</b> | <b>90sec</b> |                          |
| <b>Urokinase</b>                      | <b>20uL</b>  |             | <b>60sec</b> |                          |
| <b>Chromogenic substrate in HN</b>    | <b>100uL</b> |             | <b>10min</b> | <b>Color development</b> |
| <b>Acetic acid</b>                    | <b>100uL</b> |             |              | <b>Stop reaction</b>     |

4. Plate is read at 405 nm
5. Absorbance are converted to %NMP plasminogen using the standard curve