

Protein C Assay

Measured by the level of chromogenic substrate cleaved by activated protein C

Materials

1. Versa Max microtiter plate reader (Molecular Devices, CA)
2. 96-well microtiter plate
3. Citrated plasma samples
4. Normal mouse plasma (NMP)
5. TBS with 100 mM CsCl
6. 2 U/mL Protein C activator
7. 2.5 mM chromogenic substrate specific for activated Protein C (APC)
8. 20% acetic acid

Procedure

1. Plasma samples may be diluted to 1:10 - 1:20 in TBS-CsCl
2. Standard curve is prepared with each plate by diluting NMP 1:4 to 1:64 in TBS, analyzed simultaneously on the sample plate
3. Plate wells are loaded as follows:

Sample dilutions, in duplicate	10uL	37°C	15min	
Protein C activator	25uL	37°C	15min	
Chromogenic substrate	25uL	37°C	1hr	Color development; plate must be covered
Acetic acid	25uL	37°C	1hr	Stop reaction

4. Plate is read at 405 nm
5. Absorbance are converted to %NMP Protein C using the standard curve