

Protein S Assay

Protein S antigen is detected by binding to an antibody, then secondary binding by a conjugated antibody that will produce a detectable color change with the addition of a substrate

Materials

1. Versa Max microtiter plate reader (Molecular Devices, CA)
2. 96-well microtiter plate
3. Citrated plasma samples
4. Normal mouse plasma (NMP)
5. 10 µg/µL rabbit anti-human Protein S polyclonal antibody in 50 mM Na₂CO₃, pH 9.6
6. TBS
7. 3% BSA in TBS
8. 1% BSA in TBS
9. 0.05% Tween 20 in TBS
10. Horseradish peroxidase (HRP)-conjugated rabbit anti-human Protein S polyclonal antibody diluted 1:1000 in 1% BSA in TBS
11. 1 N H₂SO₄

Procedure

1. Plasma samples may be diluted to 1:100 - 1:200 in TBS/1% BSA
2. Log-log standard curve is prepared with each plate by diluting NMP 1:25 to 1:100 in TBS/1% BSA, analyzed simultaneously on the sample plate
3. Plate wells are loaded as follows:

Rabbit anti-human protein S Ab	100uL	5°C	Overnight	
3%BSA in TBS	200uL	37°C or 5°C	3-5hr or Overnight	Block
1%BSA in TBS				Wash 1x
Sample dilutions, in duplicate	100uL	5°C	Overnight	
0.05% Tween 20 in TBS				Wash 5x
Diluted HRP-conjugated rabbit Ab in TBS/1%BSA	100uL	5°C	Overnight	
0.05% Tween 20 in TBS				Wash 5x
TMB peroxidase substrate			4hr	Color development
H2SO4	100uL			Stop reaction

4. Plate is read at 655 nm during development and then at 450 nm after reaction is stopped
5. Absorbance are converted to %NMP Protein S using the standard curve